



---

## Quick Reference Guide





## Opus Plasma : GLIDE Tip Treatment Guidelines

### Recommended Treatment Parameters

Skin Type*	Energy (W)	Plasma Intensity (%)	Passes (#)
I – III	15 – 20	40	4
IV – VI	12 – 17	40	4

### Pre-Treatment Guideline:

- Apply topical anesthetic prior to treatment.
- Thoroughly clean and dry treatment area.
- Remove all jewelry from treatment area.
- Always perform a skin test on the intended treatment area (Table 1).
- Assess skin condition after 5 minutes.

Mild to moderate erythema should appear.

Do not begin the treatment if other skin responses appear (e.g. irritation, allergic reaction).

Skin Type*	Energy (W)	Plasma Intensity (%)	Passes (#)
I – III	15 – 20	40	1
IV – VI	12 – 17	40	1

Table 1. Skin Test Parameters for Glide Tip

### Treatment Technique and Setting Guideline:

- Set recommended treatment parameters and press “standby” on screen.
- Position the handpiece perpendicular to the skin first before pressing the foot-switch/trigger to avoid premature energy discharge when tip comes in close contact to the skin.
- Treatment is carried out by several passes over the treatment area in a one-directional, crosshatched (X-Y) pattern.
- Keep the speed of treatment consistent to 3 – 4 cm/sec.  
A consistent, gentle buzz is an audible guide to ensure consistent Plasma delivery.
- 4 passes are recommended for about 10% ablation coverage.
- Endpoint: Mild to moderate erythema and edema. Superficial epidermal ablation zones may be visible.
- 2 – 3 treatments are recommended.
- Intervals: 3 – 4 weeks. Treat only on healed skin.
- Glide tip are single-use only.

### Post-Treatment Guideline:

- Cool the skin with cold packs or cold air.
- Clean treated area, use a skin calming product if necessary and follow post-operative care guidelines.
- Instruct patient to avoid hot showers 24 hours post procedure.
- Avoid sun exposure after and in between treatments and use UVA/UVB sunscreen (SPF>30).



## Opus Plasma : FOCUS:IN / FOCUS Tips Treatment Guidelines

### Recommended Treatment Parameters

Skin Type*	Energy (W)	Pulse Duration (sec)	Plasma Intensity (%)	Passes (#)
I – III	40 – 50	1.0	25	1 – 2
IV – VI	35 – 40	1.0	25	1 – 2

### Pre-Treatment Guideline:

- Apply topical anesthetic prior to treatment.
- Thoroughly clean and dry treatment area.
- Remove all jewelry from treatment area.
- Always perform a skin test on the intended treatment area (Table 2).

Skin Type*	Energy (W)	Pulse Duration (sec)	Plasma Intensity (%)	Passes (#)
I – III	40 – 50	0.5	25	1
IV – VI	35 – 40	0.5	25	1

Table 2. Skin Test Parameters for FOCUS:IN / FOCUS Tips

- Assess skin condition after 5 minutes. Mild to moderate erythema should appear. Do not begin the treatment if other skin responses appear (e.g. irritation, allergic reaction).

### Treatment Technique and Setting Guideline:

- Set recommended treatment parameters and press “standby” on screen.
- Position the FOCUS:IN / FOCUS tip gently on the skin and press the foot-switch/trigger to deliver energy. Move tip to the next spot and press trigger again.
- Treatment may be carried out by more than one pass over the treatment area.
- Endpoint: Mild to moderate erythema and edema.
- 2 – 3 treatments are recommended.
- Intervals: 3 – 4 weeks. Treat only on healed skin.
- FOCUS:IN and FOCUS tips are reusable and cleaning is recommended before the next session. Refer to the Operator Manual for recommended cleaning instructions.

### Post-Treatment Guideline:

- Cool the skin with cold packs or cold air.
- Clean treated area, use a skin calming product if necessary and follow post-operative care guidelines.
- Instruct patient to avoid hot showers 24 hours post procedure.
- Avoid sun exposure after and in between treatments and use UVA/UVB sunscreen (SPF>30).